

USE OF CYANOACRYLATE ADHESIVES

To make teeth additions or relines in Deflex polymers that do not have chemical adhesion to acrylic, the use of cyanoacrylate-based adhesives is recommended as a bonding agent that contributes to adhesion with acrylic.

Prior to using a cyanoacrylate-based adhesive for the reparation of a dental denture, the use information provided by the adhesive manufacturer should be read to verify that it is suitable for permanent contact with the oral mucosa.



Technical department

info@deflex.com.ar - www.deflex.com.ar