

# Flyer BulkFlex FRENTE 14 cm x 21 cm



## Injectable materials in bulk

Aesthetic, light-weight and resistant dentures

**BULKFLEX** is Deflex's new line of injectable materials in bulk.

High quality thermoplastic polymers in bulk packaging allows for dispensing only the amount necessary for each denture.

- Resistant materials
- Metal-free
- Minimal shrinkage
- High color stability
- Easy injection



### **TENAZ** Semi-flexible polyamide

- Highly resistant to breaking or deformation
- Very fluid, easy injection
- Translucency: medium

*Recommended for flexible partials*



### **VERSATIL** Semi-rigid polyamide

- Total versatility: excellent balance of elasticity and strength
- Allows for all types of design. Occlusal rests can be made
- High translucency. Easy to polish

*Ideal for all types of partials*



### **ACRYNEL** Injectable acrylic

- Extra strength translucent acrylic
- Minimal shrinkage
- Adheres to normal acrylic (possibility of relining and additions using conventional materials and techniques)

*Light and resistant full dentures*



### **ACETAL** Injectable Acetal

*Partial denture frameworks and tooth-colored retainers*

#### Uses:

- Partial dentures
- Full dentures
- Combined with chromium
- Overdentures
- Occlusal splints
- Partial denture frameworks
- Tooth-colored retainers

**NUXEN**  
SRL

 Ayacucho 1053; 3°A - City of Buenos aires - Argentine

 (+5411) 4812-9638

 (+54911) 11-6720-5702

 info@deflex.com.ar

ET 62 Rev00 01-01-2017